

HISTO/Orientator™



“De-wrinkle” your paraffin sections



“De-wrinkle” your paraffin sections with TBS’

HISTO/Orientator™



The convenient 2” x 2” heated surface can be used for both “**de-wrinkling**” paraffin sections or for the “**rapid attachment**” of fresh frozen specimens.

A quick touch of the slide to the hot plate will remove wrinkles and bubbles on your sections. You can also keep a unit on your cryostat to ensure solid attachment of frozen sections to slides to prevent loosening of the tissue during staining agitation.

Ordering Information:

Cat #	U/M	Description
TO-120	ea	HISTO/Orientator™, w/ 2”x2” heated surface; 100-120 VAC, 50/60Hz, 0.5A; 2”W x 2”D x 3”H
TO-220	ea	HISTO/Orientator™, w/ 2”x2” heated surface; 220-240 VAC, 50/60Hz, 0.25A; 2”W x 2”D x 3”H



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com

Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability. HISTO/Orientator is a trademark and TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. © February 2006

Bulletin # 11.1.15