

Flotation Workstation



Two Smart Designs from TBS[®]


TBS[®]

Equip your lab with a TBS Flotation Workstation

and make your tissue preparation easier, safer and more efficient.

Imagine...


the practicality of a removable glass dish, the long term aesthetics afforded by a scratch and chemical resistant plastic housing, the optimum specimen viewing provided by a high contrast background and the dependability & safety of microprocessor controlled electronics ... **combined with the convenience of a HISTO/Orientator™ for flattening sections and a slide dryer all in one compact unit.**



For those flotation bath applications not requiring the added convenience of a slide dryer or HISTO/Orientator™ the basic TBS flotation bath provides all the functional designs of the Flotation Workstation at an affordable price.

Ordering Information:

Cat #	U/M	Description
FWS-120	ea	Flotation Workstation, microprocessor controlled w/ removable glass bowl, 2" x 5" convection slide dryer and 2" x 2" HISTO/Orientator™; 100-120 VAC, 50/60Hz, 3.3A; 14"W x 18"D x 5.25"H; 10.35lbs
FWS-220	ea	Flotation Workstation, microprocessor controlled w/ removable glass bowl, 2" x 5" convection slide dryer and 2" x 2" HISTO/Orientator™; 220-240 VAC, 50/60Hz, 1.75A; 14"W x 18"D x 5.25"H; 10.35lbs
TFB-120	ea	Tissue Flotation Bath, microprocessor controlled w/ removable glass bowl and built-in clock; 100-120 VAC, 50/60Hz, 3.3A; 14"W x 18"D x 5.25"H; 10lbs
TFB-220	ea	Tissue Flotation Bath, microprocessor controlled w/ removable glass bowl and built-in clock; 220-240 VAC, 50/60Hz, 1.75A; 14"W x 18"D x 5.25"H; 10lbs
RB	ea	Replacement Bowl, 8(20) x 8(20) x 2(5) in(cm)
RFL-9W	pk	Fluorescent lamp, replacement for FWS and TFB; 2ea/pk



Please visit our website for additional product information, technical bulletins, career opportunities, MSDS files, etc.

www.trianglebiomedical.com

Distributed By:

TBS TRIANGLE BIOMEDICAL SCIENCES, INC.

3014 CROASDAILE DRIVE • DURHAM, NORTH CAROLINA 27705-2507 • USA • 919.384.9393 • FAX:919.384.9595 • EMAIL: tbs@trianglebiomedical.com
 Specifications are current at time of publication; however, they are subject to change without notification. Please contact our North Carolina office or one of our distributors for current prices and availability.
 HISTO/Orientator is a trademark and TBS is a registered trademark of Triangle Biomedical Sciences, Inc., Durham, North Carolina, USA.

Litho in U.S.A. © April 2006

Bulletin # 8,9,2.0